

Transesophageal Electronic Sector Probe

UST-52126

Instruction Manual

MN1-5791 Rev.6

*Notes for operators and responsible maintenance personnel*

- ★ *Please read through this Instruction Manual carefully prior to use.*
- ★ *Keep this Instruction Manual together with the ultrasound diagnostic instrument for any future reference.*



 **Hitachi, Ltd.**

© Hitachi, Ltd. 2013, 2017. All rights reserved.



## Introduction

This is an instruction for model UST-52126, an ultrasound probe.

Read the manual carefully before using the instrument. Take special note of the items in section 1, "Safety Precautions".

Keep this manual securely for future reference.

The CE mark on the probe indicates that this probe is valid when it is connected to equipment bearing the CE mark that is specified as available in section 2 of this document. Therefore, if a probe bearing the CE mark is connected to equipment that is specified as available but does not have a CE mark, part of this instruction manual may not apply.

## Symbols used in this document

The terms below are used in the safety information provided to prevent hazards and injuries to the operator or patients. The severities of the hazard and injury that can occur when failing to observe the displayed safety information are indicated in four levels: "Danger", "Warning", "Caution" and "Note".

 <b>Danger</b>
-------------------------------------------------------------------------------------------------

Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury to the operator or patient.
---------------------------------------------------------------------------------------------------------------------------------------

 <b>Warning</b>
----------------------------------------------------------------------------------------------------

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury to the operator or patient.
----------------------------------------------------------------------------------------------------------------------------------------

 <b>Caution</b>
----------------------------------------------------------------------------------------------------

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the operator or patient, or property damage only.
----------------------------------------------------------------------------------------------------------------------------------------------------------------

 <b>Note</b>
-------------------------------------------------------------------------------------------------

Indicates a strong request concerning an item that must be observed in order to prevent damage or deterioration of the equipment and also to ensure that it is used efficiently.
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

The type of safety information is indicated by the symbols below.

 This symbol means attention is required.
------------------------------------------------------------------------------------------------------------------------------

 This symbol means that the described action is prohibited.
------------------------------------------------------------------------------------------------------------------------------------------------

 This symbol means the described action is mandatory.
------------------------------------------------------------------------------------------------------------------------------------------

# CONTENTS

1. Safety Precautions .....	1
1-1. Intended use .....	1
1-2. Usage precautions.....	1
1-2-1. Warnings and safety information .....	2
1-2-2. Transducer cover usage precautions .....	5
1-2-3. Cleaning, disinfection and sterilization precautions .....	6
1-2-4. Labels.....	7
2. Specifications and Parts name.....	11
2-1. Principles of operation .....	11
2-2. Specifications .....	12
2-2-1. Specifications of the probe.....	12
2-2-2. Specifications of the bite block.....	13
2-2-3. Transducer cover .....	13
2-3. Performance .....	14
2-4. Name of each parts .....	14
2-5. Environmental conditions .....	16
2-5-1. Operating environmental conditions .....	16
2-5-2. Storage environmental conditions.....	16
2-6. Classification of ME equipment .....	16
3. Preparations for Use .....	17
3-1. Start up check .....	17
3-1-1. Visual check .....	17
3-1-2. Deflection portion operation check .....	17
3-1-3. Transducer rotation part operation check.....	17
3-1-4. Verification of cleaning, disinfection and sterilization.....	18
3-1-5. Verification of operation .....	18
4. Usage.....	19
4-1. Operation.....	21
4-1-1. Operation of the bite block.....	21
4-1-2. Transducer cover usage precautions.....	22
4-1-3. Removal of transducer cover .....	23
4-1-4. Manipulation of the deflection portion.....	24
4-1-5. Manipulation of the transducer rotation switch .....	25
4-1-6. Manipulation of the Lock knob and Lock lever.....	26
4-1-7. Temperature Control System .....	27
4-1-8. Removing the insertion portion .....	27
4-2. Connecting to the ultrasound diagnostic instrument .....	28
4-3. Removing from the ultrasound diagnostic instrument.....	29

4-4. Actions to be taken when an abnormal state is detected.....	30
4-4-1. Ensuring safety of patients.....	30
4-4-2. Handling the instrument.....	30
5. Cleaning, disinfection and sterilization .....	31
5-1. Precautions for cleaning, disinfection and sterilization.....	32
5-2. Reprocessing instruction according to ISO 17664 .....	33
5-3. Point of use (Pre-cleaning).....	35
5-4. Containment and transportation.....	35
5-5. Manual cleaning and disinfection .....	36
5-5-1. Manual cleaning .....	37
5-5-2. Manual disinfection .....	38
5-5-3. Cable and connector .....	39
5-6. Automated cleaning and disinfecting .....	39
5-7. Applicable cleaners and disinfectants / Suppliers List .....	40
5-8. Drying.....	42
5-9. Maintenance, inspection and testing.....	42
5-10. Packaging .....	42
5-11. Sterilization.....	43
5-11-1. Ethylene oxide (EtO) gas sterilization .....	43
5-11-2. STERRAD® sterilization .....	44
5-11-3. Liquid sterilization (USA only) .....	45
5-11-4. Autoclave sterilization.....	45
5-11-5. STERIS® sterilization.....	46
5-12. Storage.....	46
6. Storage .....	47
6-1. Actions before storing the probe.....	47
6-2. Environmental conditions for storage.....	47
7. Moving and Transporting.....	49
7-1. Moving and transporting.....	49
7-2. Preparing the probe and accessories for moving .....	49
7-3. Packing for transportation.....	49
7-4. Environmental conditions during transportation .....	49
8. Periodic Inspection.....	51
8-1. Safety tests.....	51
8-2. Testing of measurement tolerances.....	52
8-2-1. Conducting tests .....	52
8-2-2. Result judgment.....	52

9. Configuration .....	53
9-1. Standard configuration .....	53
9-2. Transducer cover .....	53
10. Disposal of the Device.....	55

This Instruction Manual contains the main body of 56 pages and 6 pages until the CONTENTS.

## 1. Safety Precautions

### 1-1. Intended use

This probe is intended for use by a doctor when inserted into the patient's esophagus and its tip contacts the esophageal wall making ultrasonic observations of the heart.

Please refer to the ultrasound diagnostic instrument instruction manual used with this probe for the probe intended use information.

Regarding with the connectable instrument, please refer to section 2-2. Specifications of this manual.

#### Caution

 Do not use this equipment for other than its intended purpose.  
Use for other purposes can cause burns or other injuries to the patient or operator.

### 1-2. Usage precautions

The terms below are used in the safety information provided to prevent hazards and injuries to the operator or patients. The severities of the hazard and injury that can occur when failing to observe the displayed safety information are indicated in four levels: "Danger", "Warning", "Caution" and "Note".

#### Danger

Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury to the operator or patient.

#### Warning

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury to the operator or patient.

#### Caution

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the operator or patient, or property damage only.

#### Note

Indicates a strong request concerning an item that must be observed in order to prevent damage or deterioration of the equipment and also to ensure that it is used efficiently.

The type of safety information is indicated by the symbols below.

 This symbol means attention is required.

 This symbol means that the described action is prohibited.

 This symbol means the described action is mandatory.


1-2-1. Warnings and safety information

 <b>Warning</b>	
	Operation must be performed by a skilled doctor. Improper operation can injure the patient. Operations must be performed by a doctor who fully understands the characteristics of ultrasound diagnostics and who is skilled and has a thorough knowledge of transesophageal echocardiography.
	Do not use for patients who have esophageal barriers. This may cause serious injury to the patient.
	Follow the information in this manual and the documentation supplied with any equipment used together with this probe. Use that is not in accordance with the supplied documentation can result in a serious or moderate injury, equipment breakdown, or physical damage that impairs operation.
	Be sure to preparations for use. Use of the equipment while failing to notice an abnormal condition can result in injury to the operator or patient. If any abnormalities are noted on the equipment in the start up check, immediately stop using it and contact one of our offices and/or distributor's offices listed on the back cover. See section 3-1 "Start up check" for the start up check content and procedure.
	Do not use on the eyes. This probe is not intended for use on the eyes. The acoustic output can have an adverse effect on the eyes.
	Do not attempt to disassemble, modify, or repair the equipment. Electric shock or other unforeseen accidents could result. Contact one of our offices and/or distributor's offices listed on the back cover to request repair.
	Do not immerse the probe connector and proximal head into water or chemical solution, it will cause a breakdown of the probe. The probe is not waterproof. Regarding information about the waterproof part of the probe refer to chapter 1-2-3 and chapter 5.
	Clean, disinfect and sterilize before using the probe. Perform proper cleaning, disinfection and sterilization after use. Otherwise, there is a risk of infection. Note that the probe is not sterilized at the factory. Before using the probe first, be sure to clean, disinfect and sterilize it.
	Wear medical gloves during examination. Conducting examinations with the bare hands can expose the operator to a risk of infection.
	This probe must not be used in direct with the heart. This may cause patient to receive an electric shock.
	Select the size of the probe with the patient's physical constitution in mind. Using a portion that is too large and unfit for the patient's physical constitution is very dangerous and can harm the patient.
	Use the transducer cover over the insertion portion. If the transducer cover is not used and will contact Lignocaine Hydrochloride such as xylocaine jelly or gel , this may cause exterior deterioration. If the transducer cover is not used, residual pathogens on the probe could infect the patient.
	Use the sterilized echo jelly (gel) which is attached to the transducer cover. When using an acoustic medium which is not recommended in the instruction manual, it causes deterioration of the probe surface.
	Use the bite block supplied with the probe. If the probe is bitten, the probe may be damaged and a hazard to the patient will occur.


**⚠ Warning**

**!** Do not try to forcibly perform operations.  
Excessive force causes injury to the patient. If an abnormal resistance force is felt, stop use of the equipment.


**!** If the probe gets bent in the esophageal (J turn), do not try to move the probe by force. If the image reverses suddenly or the image disappears suddenly at the time of the insertion, J-turn should occur. If a J-turn should occur, do not try to pull back the probe by force, but rather, preferably under x-ray monitoring to check the state of the probe, carefully push the probe until its tip comes into the stomach and straighten the bend by operating the angle knob before pulling back the probe. Forcing it to move can injure the patient.



Occuring J-turn



The probe gets bent  
in the esophageal



Push the probe until stomach  
and straighten

**!** Keep the angle knob in the free position during an operation of pulling out.  
Removing the deflection portion while locked may injure the patient.

**!** Dispose the probe used for patients with Creutzfeldt-Jakob disease.  
Otherwise, there is a risk of infection to the operator or patient. Our ultrasound probe is not compatible with any disinfection/sterilization method for Creutzfeldt-Jakob disease.

**⊘** Do not use the equipment fallen on to floor.  
Otherwise, there is a risk of infection. Stop the operation and perform the procedure in section 8 "Periodic Inspection", section 5 "Cleaning, disinfection and sterilization" and section 3-1 "Start up check".

 <b>Caution</b>	
	Constantly check for anything abnormal about the patient's condition and the equipment. Continued use without noticing that an abnormal condition has occurred can result in an electric shock and injury to the operator or patient. If an abnormal condition occurs, immediately move the equipment away from the patient and stop use of the probe.
	The equipment is vulnerable to damage by impact. Therefore, handle it with care. There is a risk of damage to the equipment when the equipment is fallen or hit somewhere.
	Do not use this probe with other equipment except for those specifically approved in the manual. Use with unapproved equipment can result in an electric shock, burn, or other injury to the patient or operator and damage to the probe and the other equipment.
	Scan for the minimum length of time necessary for the diagnosis and at the lowest suitable output. Overuse can adversely affect the internal tissues of the patient. For details about the acoustic output, please refer to the documentation supplied with the ultrasound diagnostic instrument.
	Constantly check about surface temperature of probe tip. Do not exceed the maximum surface temperature of probe tip by acoustic power. See section 4-1-7.
	Regularly perform maintenance inspection and safety tests of the equipment. If you use equipment for a long period of time, it can reduce the performance, or cause smoke or fire. If anything unusual occurs, immediately stop using it and contact one of our offices and/or distributor's offices listed on the back cover.
	Use, move and transport the equipment under the environmental conditions specified in this manual. Otherwise, it may be damaged. See section 2-5 "Environmental conditions" and section 7-4 "Environmental conditions during transportation".
	Do not use with a defibrillator. This could result in the probe to be damaged. If you have to use a defibrillator with the probe in place, make sure to conduct a safety tests according to Section 8.
	Do not bend or twist the deflection portion unnecessarily or manually. This could make the probe unusable.
	Do not bend or pull the insertion tube unnecessarily. This could make the probe unusable.
	Keep the angle knob in the free position during an operation of curving. Operating the probe while locked, it may become unusable.
	Do not apply olive oil or xylocaine spray , xylocaine (Lignocaine Hydrochloride) jelly or gel to the probe. This may cause exterior deterioration.

## 1-2-2. Transducer cover usage precautions

 <b>Warning</b>	
	Use the transducer cover over the insertion portion. If the transducer cover is not used and will contact Lignocaine Hydrochloride such as xylocaine jelly or gel , this may cause exterior deterioration.
	Use Hitachi-approved transducer cover only. If a transducer cover is used which is not recommended by Hitachi, it can cause an adverse reaction to the body of the patient.
	Check that the transducer cover is sterilized. Use of an infective item could spread infection to the patient.
	Use the sterilized echo jelly (gel) inside the transducer cover which is attached to the transducer cover. When using an acoustic medium which is not recommended in the instruction manual, it causes deterioration of the probe surface.
	Do not reuse the transducer cover . Use of an infective item could spread infection to the patient.
	Do not apply unsterilized acoustic medium to the outer surface of the transducer cover . Use of an acoustic medium that is contaminated by a pathogen can cause an infection on the patient.
	Do not use on patients who may have an allergic reaction to latex products. Use of the transducer cover for these types of patients could result in anaphylactic shock. Ask the patient about allergy history beforehand.


 <b>Caution</b>	
	Check the transducer cover for abnormalities before use. Regarding the storage of the transducer cover, follow the instructions on the package of the transducer cover.
	Check that the acoustic medium has no air bubbles inside the transducer cover that is covering the probe. Air bubbles inside the transducer cover can result in misdiagnosis caused by overlooking or misinterpreting lesions due to poor image quality or improper rendering.

1-2-3. Cleaning, disinfection and sterilization precautions


 <b>Warning</b>	
	Wear protective gloves and other protective gear during cleaning, disinfection and sterilization. Handling of the equipment with your bare hands before disinfection or sterilization can result in an infection.
	After soaking in cleaning agents, thoroughly wash the equipment with running water. Residual cleaning agents can cause an adverse reaction on the bodies of the operator or patient.
	After soaking in a disinfectant, thoroughly wash the equipment with deionized water. Residuals of the disinfectant can cause an adverse reaction on the bodies of the operator or patient.
	Perform aeration completely after gas disinfection and sterilization. Residual gas can cause an adverse reaction on the bodies of the operator or patient.
	Do not clean, disinfect or sterilize using procedures other than those specified in this manual. Infection could result due to incomplete cleaning, disinfection and sterilization. It can also result in damage to the equipment or reduced performance. The probe cannot withstand autoclave sterilization or boiling and other types of sterilization at temperatures exceeding 60°C (140°F).
	Do not perform gas disinfection and any sterilization. This probe cannot withstand Gas disinfection, and any types of sterilization.
	For details on the usage conditions of chemicals and sterilization procedures, refer to the documentation supplied with the respective chemical or sterilization equipment. Infection could result due to incomplete disinfection or sterilization. This could also cause deterioration of the equipment.

 <b>Caution</b>	
	Do not immerse the probe in any liquids beyond the range of IPX7 shown in section 2-2 "Specifications". Don't submerge or immerse Proximal head and Connector into water or chemical solution. Use when liquid has gotten inside the connector can result in a risk of electric shock to the operator or patient. If liquid gets inside the connector, immediately stop use and contact one of our offices and/or distributor's offices listed on the back cover.


1-2-4. Labels  
(1) Probe unit



Label 1




Label 2



Label 3

Label 1



Electronic sector probe mark

Frequency

Label 2



This instrument complies with Directive 93/42/EEC relating to Medical Device and Directive 2011/65/EU relating to RoHS.



IPX7 mark  
See section 2-2, "Specifications".



Type BF applied part



Do not waste the instrument as general waste. Comply with a local regulation.  
See section 10.



Safety warning sign



Biohazard  
See section 5.



Follow the instruction manual to operate this instrument.  
If not avoided, may result in injury, property damage, or the equipment trouble.

Label 3(Example)



Country of manufacture

Model

Serial No.



Manufacturer

Address


Rx Only:

By prescription only. U.S. Federal Law restricts this device to sale on order of a physician only.

(2) Storage case



Label A



Model  
Serial No.

Label B



This instrument complies with Directive 93/42/EEC relating to Medical Device and Directive 2011/65/EU relating to RoHS.



DATE OF MANUFACTURE  
(in case of 2016-09)




MANUFACTURER



AUTHORISED REPRESENTATIVE IN EUROPEAN  
COMMUNITY

Label C



DATE OF MANUFACTURE  
(in case of 2016-09)

---

## 2. Specifications and Parts name

---

### 2-1. Principles of operation


This probe and the ultrasound diagnostic instrument enable image diagnosis using ultrasonic waves. These instruments operate under the principles described below.

- (1) When an electric pulse signal is applied from the transmitter to the transducer of the probe, the transducer operates by converting electrical vibrations to mechanical vibration energy for emitting pulse-shaped ultrasonic waves into the body part contacting the transducer or into liquid or other acoustic medium.
- (2) The emitted ultrasonic waves are reflected by boundaries with different acoustic characteristics (acoustic impedance) within the body.
- (3) The transducer is also used to receive reflected ultrasonic waves. The transducer vibrates mechanically due to the received ultrasonic vibrations and uses an electro-mechanical conversion operation to convert the received mechanical vibrations to electric energy. The received echo is also converted to electric signals and a brightness modulation operation is used to convert the electric pulses to shades of brightness for forming an image.

## 2-2. Specifications


### 2-2-1. Specifications of the probe

Application regions:	Heart and thoracic aorta
Form of application to patient:	Transesophageal
Connectable instruments:	Prosound F75, ARIETTA 70, ARIETTA 60
	NOTE:
	At the time of publication of this manual, the connectable diagnostic ultrasound instrument or instrument software version available with this probe is different for each country, please refer to the instrument instruction manual or contact your local Hitachi representative.
Field of view:	90°
Scan direction:	Any direction in ±90° relative to insertion direction of insertion portion
Acoustic Working Frequency:	5.1±20% MHz
Direction and range of deflection:	UP120° DOWN 90° LEFT 45°RIGHT 45°
Outer diameter of insertion portion:	φ 11mm
Effective insertion distance:	1.1m
Cable length:	1.56 m
Weight:	1500 g
Service life:	Three years
Range of applied part:	As shown in the figure below.
Parts treated as applied parts:	As shown in the figure below.
IPX7 range:	As shown in the figure below.
External dimensions:	As shown in the figure below
Maximum surface Temperature of Probe tip	43°C



2-2-2. Specifications of the bite block

Material: Polyetherimide  
Service life: Three years  
External dimensions: As shown in the figure below.



Remarks  
The dimensions are within  $\pm 10\%$  of the indicated values.

2-2-3. Transducer cover

Use by covering the transducer cover over the insertion portion.


Transducer cover	CIVCO Transducer cover	610-933
------------------	------------------------	---------

Remarks  
This transducer cover is not included in this probe-kit.

2-3. Performance

For measurement tolerances, operating tolerances and other data, refer to the instruction manual for the ultrasound diagnostic instrument.

2-4. Name of each parts



 <b>Caution</b>	
	Do not pull, bend, twist, or apply excessive force to the cable. The conductors may break and the cable may become unusable.
	Do not subject the ultrasonic irradiation area to hard impact. The impact may cause damage to the transducer, and that results in noise or no echo in the image. In most cases, the ultrasonic irradiation area itself is not damaged because the part is made of plastic material.
	Do not bend or twist the deflection portion unnecessarily or manually. This could make the probe unusable.
	Keep the angle knob in the free position during an operation of curving. Operating the probe while locked, it may become unusable.
	There is no mark to indicate front direction of probe. It becomes a right angle operation when continuing pushing the transducer rotation switch and turning the transducer completely. See section 4-1-5.

**Remarks**

According to the property of a acoustic medium which filled a tip ,the change of color may be seen in membrane.  
But there is no influence on the performance and safety of the probe.

## 2-5. Environmental conditions

Use and store the equipment under the following conditions.

### 2-5-1. Operating environmental conditions

Ambient temperature:	10°C to 40°C 50°F to 104°F
Relative humidity:	30% to 75%
Atmospheric pressure:	700 hPa to 1060 hPa
Altitude:	3,000 m or less

### 2-5-2. Storage environmental conditions

Ambient temperature:	-10°C to 50°C 14°F to 122°F
Relative humidity:	10% to 90%
Atmospheric pressure:	700 hPa to 1060 hPa

#### Caution

 Avoid operating or storing the equipment in the following locations.

- Locations exposed to water or other liquids
- Locations subject to adverse conditions such as air pressure, temperature, humidity, ventilation, direct sunlight, dust, or air containing salt, sulfur, or other corrosive substances
- Locations where chemical substances are stored or where gases are generated

Storage in these locations can result in a breakdown or reduced performance.

 Avoid rapid temperature change which may cause condensation. Avoid using in locations where condensation or water droplets can form.

Condensation can occur when moving the probe from a cool location to a warm one. Use when condensation has occurred can result in a breakdown or reduced performance.

## 2-6. Classification of ME equipment

- Classification based on degree of protection against electric shock . Type BF applied part
- Classification for protection against ingress of liquids ..... IPX7 (Watertight equipment)
- Operation mode ..... Continuous operation
- Method of sterilization ..... See section 5 “Cleaning, disinfection and sterilization”

For the range of applied parts, parts treated as applied parts and the range of IPX7, see section 2-2-1.

## 3. Preparations for Use

### Warning



Be sure to preparations for use.

Using the equipment without noticing an abnormal condition can result in injury to the operator or patient. If an inspection finds an abnormal condition in the equipment, immediately stop use and contact one of our offices and/or distributor's offices listed on the back cover.

### 3-1. Start up check

#### 3-1-1. Visual check

Visually and tactually check and confirm the following:

- Insertion portion is free from holes, dents, scratches, cracks, deformations, color changes and the like on the surface. Also there is no leakage of the oil inside the transducer.
- Manipulate the deflection portion using the angle knob and check for any protrusions or cracks in this area.
- Lightly grip insertion tube and deflection portion by hand and let it slide. Then, it shall neither catch your hand on the way nor shall it be slack.
- Make sure that all surfaces of the probe connector, cable and proximal head are free of scratches, cracks or exfoliation.
- Bite block is free from holes, dents, scratches, cracks, deformations, color changes and the like on the surface.
- Transducer cover is free from holes, dents, scratches, cracks, deformations, color changes and the like on the surface.

#### 3-1-2. Deflection portion operation check

- Gently turn the angle knob in each direction until it stops and check the following:
  - \* Make sure there are no catching or irregularities in force to turn the angle knob.
  - \* Make sure the deflection portion is bent smoothly in all directions.
- Operate the lock lever to make sure there are no abnormalities in the curvature holding or releasing functions.

#### 3-1-3. Transducer rotation part operation check

- Push the switches in each direction of the transducer as far as it will go on the displayed screens, and confirm the following.
  - \* Make sure the transducer rotates smoothly in all directions.

### 3-1-4. Verification of cleaning, disinfection and sterilization

Verify that cleaning, disinfection and sterilization are conducted according to the intended use.

### 3-1-5. Verification of operation

Connect to the ultrasound diagnostic instrument by following the instructions in section 4-2, "Connecting to the ultrasound diagnostic instrument" and check that the selected probe match the sector display and the displayed frequency and check the image for errors.

#### Remarks

For details on the displayed screens, see the documentation supplied with the ultrasound diagnostic instrument.

If the probe is operated in still air, brightness on the top of the image may be non uniform, but this does not affect the performance of the probe.

 <b>Caution</b>
 Do not use the probe if the selected probe and image do not match the frequency. An incorrect acoustic output can result in burns or other injuries to the patient. Contact one of our offices and/or distributor's offices listed on the back cover.


## 4. Usage

 Warning	
	Operation must be performed by a skilled doctor. Improper operation can injure the patient. Operations must be performed by a doctor who fully understands the characteristics of ultrasound diagnostics and who is skilled and has a thorough knowledge of transesophageal echocardiography.
	Do not use for patients who have esophageal barriers. This may cause serious injury to the patient.
	Follow the information in this manual and the documentation supplied with any equipment used together with this probe. Use that is not in accordance with the supplied documentation can result in a serious or moderate injury, equipment breakdown, or physical damage that impairs operation.
	Be sure to preparations for use. Use of the equipment while failing to notice an abnormal condition can result in injury to the operator or patient. If any abnormalities are noted on the equipment in the start up check, immediately stop using it and contact one of our offices and/or distributor's offices listed on the back cover. See section 3-1 "Start up check" for the start up check content and procedure.
	Do not use on the eyes. This probe is not intended for use on the eyes. The acoustic output can have an adverse effect on the eyes.
	Do not attempt to disassemble, modify, or repair the equipment. Electric shock or other unforeseen accidents could result. Contact one of our offices and/or distributor's offices listed on the back cover to request repair.
	Do not immerse the probe connector and proximal head into water or chemical solution, it will cause a breakdown of the probe. The probe is not waterproof. Regarding information about the waterproof part of the probe refer to chapter 1-2-3 and chapter 5.
	Clean, disinfect and sterilize before using the probe. Perform proper cleaning, disinfection and sterilization after use. Otherwise, there is a risk of infection. Note that the probe is not sterilized at the factory. Before using the probe first, be sure to clean, disinfect and sterilize it.
	Wear medical gloves during examination. Conducting examinations with the bare hands can expose the operator to a risk of infection.
	This probe must not be used in direct with the heart. This may causes patient to receive an electric shock.
	Select the size of the probe with the patient's physical constitution in mind. Using a portion that is too large and unfit for the patient's physical constitution is very dangerous and can harm the patient.
	Use the transducer cover over the insertion portion. If the transducer cover is not used and will contact Lignocaine Hydrochloride such as xylocaine jelly or gel , this may cause exterior deterioration. If the transducer cover is not used, residual pathogens on the probe could infect the patient.
	Use the sterilized echo jelly (gel) inside the transducer cover which is attached to the transducer cover. When using an acoustic medium which is not recommended in the instruction manual, it causes deterioration of the probe surface.
	Use the bite block supplied with the probe. If the probe is bitten, the probe may be damaged and a hazard to the patient will occur.


## ⚠ Warning

❗ Do not try to forcibly perform operations.  
Excessive force causes injury to the patient. If an abnormal resistance force is felt, stop use of the equipment.


❗ If the probe gets bent in the esophageal (J turn), do not try to move the probe by force. If the image reverses suddenly or the image disappears suddenly at the time of the insertion, J-turn should occur. If a J-turn should occur, do not try to pull back the probe by force, but rather, preferably under x-ray monitoring to check the state of the probe, carefully push the probe until its tip comes into the stomach and straighten the bend by operating the angle knob before pulling back the probe. Forcing it to move can injure the patient.



Occuring J-turn



The probe gets bent  
in the esophageal



Push the probe until stomach  
and straighten

❗ Keep the angle knob in the free position during an operation of pulling out.  
Removing the deflection portion while locked may injure the patient.

❗ Dispose the probe used for patients with Creutzfeldt-Jakob disease.  
Otherwise, there is a risk of infection to the operator or patient. Our ultrasound probe is not compatible with any disinfection/sterilization method for Creutzfeldt-Jakob disease.


⊘ Do not use the equipment fallen on to floor.  
Otherwise, there is a risk of infection. Stop the operation and perform the procedure in section 8 "Periodic Inspection", section 5 "Cleaning, disinfection and sterilization" and section 3-1 "Start up check".

⚠ Caution	
❗	Constantly check for anything abnormal about the patient's condition and the equipment. Continued use without noticing that an abnormal condition has occurred can result in an electric shock and injury to the operator or patient. If an abnormal condition occurs, immediately move the equipment away from the patient and stop use of the probe.
❗	The equipment is vulnerable to damage by impact. Therefore, handle it with care. There is a risk of damage to the equipment when the equipment is fallen or hit somewhere.
⊘	Do not use this probe with other equipment except for those specifically approved in the manual. Use with unapproved equipment can result in an electric shock, burn, or other injury to the patient or operator and damage to the probe and the other equipment.
❗	Scan for the minimum length of time necessary for the diagnosis and at the lowest suitable output. Overuse can adversely affect the internal tissues of the patient. For details about the acoustic output, please refer to the documentation supplied with the ultrasound diagnostic instrument.
❗	Constantly check about surface temperature of probe tip. Do not exceed the maximum surface temperature of probe tip by acoustic power. See section 4-1-7.
❗	Regularly perform maintenance inspection and safety tests of the equipment. If you use equipment for a long period of time, it can reduce the performance, or cause smoke or fire. If anything unusual occurs, immediately stop using it and contact one of our offices and/or distributor's offices listed on the back cover.
❗	Use, move and transport the equipment under the environmental conditions specified in this manual. Otherwise, it may be damaged. See section 2-5 "Environmental conditions" and section 7-4 "Environmental conditions during transportation".
⊘	Do not use with a defibrillator. This could result in the probe to be damaged. If you have to use a defibrillator with the probe in place, make sure to conduct a safety tests according to Section 8.
⊘	Do not bend or twist the deflection portion unnecessarily or manually. This could make the probe unusable.
⊘	Do not bend or pull the insertion tube unnecessarily. This could make the probe unusable.
❗	Keep the angle knob in the free position during an operation of curving. Operating the probe while locked, it may become unusable.
⊘	Do not apply olive oil or xylocaine spray, xylocaine (Lignocaine Hydrochloride) jelly or gel to the probe. This may cause exterior deterioration.

## 4-1. Operation


### 4-1-1. Operation of the bite block

Put this in the patient's mouth to keep the path of the probe inserted into the body.



#### 4-1-2. Transducer cover usage precautions

Use the sterilized echo jelly (gel) inside the transducer cover which is attached to the transducer cover.



⚠ Warning	
❗	Use the transducer cover over the insertion portion. If the transducer cover is not used and will contact Lignocaine Hydrochloride such as xylocaine jelly or gel, this may cause exterior deterioration.
❗	When using a transducer cover over the insertion portion, please use a transducer cover which is recommended by Hitachi. If a transducer cover is used which is not recommended by Hitachi, it can cause an adverse reaction to the body of the patient.
❗	Check that the transducer cover is sterilized. Use of an infective item could spread infection to the patient.
❗	Use the sterilized echo jelly (gel) which is attached to the transducer cover. When using an acoustic medium which is not recommended in the instruction manual, it causes deterioration of the probe surface.
⊘	Do not reuse the transducer cover. Use of an infective item could spread infection to the patient.
⊘	Do not apply unsterilized acoustic medium to the outer surface of the transducer cover. Use of an acoustic medium that is contaminated by a pathogen can cause an infection on the patient.

⚠ Caution	
❗	Check the transducer cover for abnormalities before use. Regarding the storage of the transducer cover, follow the instructions on the package of the transducer cover.
❗	Check that the acoustic medium has no air bubbles inside the transducer cover that is covering the probe. Air bubbles inside the transducer cover can result in misdiagnosis caused by overlooking or misinterpreting lesions due to poor image quality or improper rendering.

#### 4-1-3. Removal of transducer cover


1. Wrap the transducer cover in tissue paper and remove it from the probe.
2. Dispose used tissue paper and transducer cover using infection prevention procedures based on the rules of your facility.

#### Caution

-  Before disposing the equipment, disinfect or take other infection-prevention measures. Disposal of the equipment without taking the proper preventative measures can lead to infection.
-  After removing the equipment from the patient, check for anything abnormal about the transducer cover.  
If the transducer cover stays inside of the patient's body, the transducer cover can cause injury to the patient. Remove the transducer cover including probe with care. When the transducer cover stays inside of the patient's body, perform the required medical treatment.


#### 4-1-4. Manipulation of the deflection portion

After inserting the insertion portion into the patient's esophagus, flex the tube as indicated below in order to observe the target region.



4-1-5. Manipulation of the transducer rotation switch


After inserting the insertion portion of the probe into the esophagus of the patient, rotate the transducer using the method shown in the figure below, in order to observe the region of interest. The angle indication is a reference level.



① CCW rotation switch

This switch rotates the transducer in the counterclockwise (CCW) direction.

In continuing pushing this switch, transducer angle becomes 180 degrees.




- :Tranceducer front side  
The front mark (direction mark) side on the image display

② CW rotation switch

This switch rotates the transducer in the clockwise (CW) direction.

In continuing pushing this switch, transducer angle becomes 0 degrees.




- :Tranceducer front side  
The front mark (direction mark) side on the image display

#### 4-1-6. Manipulation of the Lock knob and Lock lever.

After inserting the insertion portion into the patient's esophagus, flex the portion as indicated below in order to observe the target region.


Lock of angle knob(for L/R operation)

Turn the lock knob to fix an angle of the deflection portion and opposite direction to free.



Lock of angle knob(for U/D operation)

Push the lock knob to fix an angle of the deflection portion and opposite direction to free.



#### 4-1-7. Temperature Control System

The tip of the probe contains a temperature sensor. This sensor monitors the surface temperature at the tip of the probe in order to prevent damage to the esophagus tissue.

When the surface temperature exceed limit, the message appears on the equipment monitor.

1) When the surface temperature becomes 41°C

When the surface temperature becomes 41°C, the message "The temperature is higher than 41.0°C!" appears on the equipment monitor.

2) When the surface temperature rise continued from 41°C

When the surface temperature rise continued from 41°C, the message "When the temperature rises any further, the Transmission is stopped." appears on the equipment monitor.

3) When the surface temperature becomes 43°C

When the surface temperature becomes 43°C, the message "TEE Thermal Limit Auto Cooling Mode in Progress" appears on the equipment monitor and the acoustic output is automatically stopped.

When surface temperature is below 41°C, This message disappears and restarts operation.

See the equipment manual for the detail.

#### 4-1-8. Removing the insertion portion

Turn free the lock lever and lock knob and straighten the deflection portion in order to remove the insertion portion from the patient' s esophagus.


Be sure to immediately clean and disinfect the probe after use.

#### 4-2. Connecting to the ultrasound diagnostic instrument

The lock lever of the connector moves over the range shown in the figure at right.

Align the ○ mark with the LOCK or RELEASE position and lock or release the electronic probe connecting socket of the diagnostic instrument ( probe connector ).

Connect the probe to the probe connector by following the procedure below.



- Connection procedure


The probe is connected when in one of the following states.

- The power switch is set to OFF.
- The image displayed on the ultrasound diagnostic instrument is frozen.

Before inserting the probe into the probe connector, check that the connector pins are not bent.

1. Turn the connector lock lever to align the ○ mark on the lever with the RELEASE position.
2. Firmly insert the connector into the probe connector.
3. Turn the lock lever clockwise by 1/4 turn until the ○ mark is aligned with the LOCK position.
4. Check that the connector is firmly inserted into the probe connector.

This completes connection of the probe.



#### ⚠ Caution




If there is resistance when trying to turn the lock lever when connecting the connector, do not forcibly try to connect it. Instead, correctly perform the steps for connecting the connector and firmly insert it into the probe connector.  
Forcibly turning the lever may damage the connector and the probe connector on the instrument.

### 4-3. Removing from the ultrasound diagnostic instrument

The lock lever of the connector moves over the range shown in the figure at right.

Align the ○ mark with the LOCK or RELEASE position and lock or release the probe connector.

Use the procedure below to remove the probe from the probe connector.



- Removal procedure

The probe is removed when in one of the following states.


- The power switch is set to OFF.
- The image displayed on the ultrasound diagnostic instrument is frozen.

1. Turn the connector lock lever to align the ○ mark on the lever with the RELEASE position.
2. Firmly grasp the connector unit and pull it out from the probe connector.

This completes the removal of the probe.

After use, perform cleaning, disinfection and sterilization of the probe by following the procedure in section 5 "Cleaning, disinfection and sterilization".

If the probe will not be used for an extended period of time, store it by following the instructions in section 6 "Storage".



#### 4-4. Actions to be taken when an abnormal state is detected

##### 4-4-1. Ensuring safety of patients

Immediately move the equipment away from the patient and quit operation.

Keep the patient in safe condition and administer the required medical treatment.

##### 4-4-2. Handling the instrument

Turn off the ultrasound diagnostic instrument, remove its plug from the AC socket and sterilize if it is contaminated. For details, refer to the instruction manual for the ultrasound diagnostic instrument.

 <b>Caution</b>
 Do not use a equipment where a problem has been found. Using a equipment in an abnormal state can cause injury to the patient. Contact one of our offices and/or distributor's offices listed on the back cover.

## 5. Cleaning, disinfection and sterilization

Applicable cleaning, disinfection and sterilization methods for each product are listed in the Table 1. The detail of each method is described in Chapter 5-2.

Table 1 Applicable cleaning, disinfection and sterilization methods

Model	Refer the corresponded items in Chapter 3, 5, 6, 8 and 11	Cleaning		Disinfection		Sterilization				
		Manual	Automated *1	Manual	Automated *1	EtO	STERRAD®	Liquid *2	Autoclave	STERIS®
UST-52126	A	X		X						
Bite block	B	X	X	X	X	X	X	X	X	X

Note: X means “Applicable”

\*1: Automated Need waterproof case (WP-001)

\*2: Liquid sterilization USA only

## 5-1. Precautions for cleaning, disinfection and sterilization

The following warnings and cautions must be observed when cleaning, disinfecting and sterilizing the probe and accessories.

 <b>Warning</b>	
	Wear protective gloves and other protective gear during cleaning, disinfection and sterilization. Handling of the probe with your bare hands before sterilization can result in an infection.
	After finishing soaking the probe in cleaning agents, thoroughly wash it with running water. Residual cleaning agents can cause an adverse reaction to the operator or the patient.
	After chemical sterilization, thoroughly wash the probe with sterile water. Residual chemicals can cause an adverse reaction to the operator or patient. (USA only)
	After disinfecting the probe, thoroughly wash the probe with deionized water. Residuals of the disinfectant can cause an adverse reaction on the bodies of the operator or patient. (EU only)
	Perform full aeration after gas sterilization. Residual gas can cause an adverse reaction to the operator or patient.
	Do not clean or sterilize using procedures other than those specified in this manual. Failure to clean and sterilize the equipment can result in an infection. It can also result in damage to the probe or reduced performance. The probe is not compatible with autoclave sterilization or boiling and other types of sterilization at temperatures exceeding 60°C [140°F].
	For details on the usage conditions of chemicals and sterilization procedures, refer to the documentation supplied with the respective chemical or sterilization equipment. Infection can be resulted due to incomplete sterilization. Wrong sterilization procedure could cause deterioration of the probe.

 <b>Caution</b>	
	Do not immerse the probe into any liquid beyond the range of IPX7. The range is indicated in the section 2-2 "specification". If any liquid enters the connector, immediately stop using the probe and contact one of our offices and/or distributor's offices listed on the back cover. Liquid in the connector could cause electric shock to the operator or patient.
	Do not wipe the ultrasonic radiation part with alcohol. Alcohol could damage the part.
	Do not use organic solvent such as thinner for cleaning to prevent the probe from damage.
	Do not use hard or sharp objects to remove residue on the probe. Such objects may damage the probe.

Additional information:

The Instructions provided above have been validated by the medical device manufacturer as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, material and personnel in the processing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

## 5-2. Reprocessing instruction according to ISO 17664

Take care about clean circumstances before using the probe on the next patients. If processors reprocess this equipment, refer to these instructions.

Table 1

WARNINGS	<ul style="list-style-type: none"> <li>• The probe is delivered unsterile. Prior to the first use, reprocess the probe.</li> <li>• Do not exceed 60 °C [140 °F].</li> <li>• Probe connector has no water resistance.</li> <li>• Particular attention is required for the cleaning of the lock lever and adjusting knob of the probe.</li> </ul>
Limitations on reprocessing	The probe is not completely submersible (Do not immerse the probe into any liquid beyond the range of IPX7. The range is indicated in the section 2-2 "specification".) Parts which are not submersible can only be disinfected by wipe disinfection.
Transportation before using	Sterile pouch or container should be kept between transportation from Central Sterile Supply Department (CSSD) to operating room. Be careful that no damages are applied to sterile pouch or container for transportation.

The level of processing required depends on the type of equipment and its use.


The CDC (Centers for Disease Control and Prevention) in the USA and the RKI (Robert Koch Institute) in Germany classify medical devices according to their use. For each classification, they specify the level of disinfection/sterilization processing that is required before use. Table 3 summarizes this information.

Table 2

Classification	Definition	Processing
Noncritical	Application part only contacts intact and uninjured skin	Cleaning ↓ Disinfection [in the USA, low-level disinfection]
Semicritical	Application part contacts mucosa (intracavitary application)	Cleaning ↓ Disinfection (Disinfectant with bactericidal, fungicidal and virucidal effect) [in the USA, high-level disinfection or sterilization]
Critical	Application part contacts intracorporeal tissue directly (intraoperative application)	Cleaning ↓ Disinfection ↓ Sterilization * <sup>1</sup>

\*<sup>1</sup>. When sterilization is not possible, the FDA in the USA recognize that disinfection (in the USA, high-level disinfection) and the use of a sterile gel and sterile transducer cover, as described in the instructions provided with the transducer cover, is an accepted method of infection control for probe.

Flowchart of reprocessing process of this probe and accessories is as follows:



NOTE: Only the accessories are compatible with automated reprocessing according to the flowchart above.

### 5-3. Point of use (Pre-cleaning)

In the operating room after use of the probe

#### A). Probe

- 1) Remove any accessories from the probe like transducer covers.
- 2) Flush patient's blood or fluid by tap water directly after use until the surface looks visually clean.
- 3) Wipe the whole surface of the probe by gauze pad and remove superficial visible impurities until the surface looks visually clean.

#### B). Bite block

- 1) Clean the bite block of all patient's blood or fluid with running tap water until the surface of the bite block looks visually clean.
- 2) Wipe the whole surface of the bite block by gauze pad and remove superficial visible impurities until the surface looks visually clean.

### 5-4. Containment and transportation

Putting the contaminated equipment into exclusive shock and damage proof container for transportation is recommended. It is recommended that instruments are reprocessed as soon as possible and not later than 4 hours after usage.

#### Caution

-  Do not immerse the probe in any liquids beyond the range of IPX7 shown in section 2-2 "Specifications". Don't submerge or immerse Proximal head and Connector into water or chemical solution. Use when liquid has gotten inside the connector can result in a risk of electric shock to the operator or patient. If liquid gets inside the connector, immediately stop use and contact one of our offices and/or distributor's offices listed on the back cover.

## 5-5. Manual cleaning and disinfection

Prepare following items before manual cleaning and disinfection.

### A). Probe

- 1) Detergent: ENZOL<sup>®</sup>/Cidezime<sup>®</sup> (Johnson & Johnson, #2258) or another cleaning agent with approved material compatibility for this medical device.
- 2) Disinfectant: Cidex<sup>®</sup> OPA (Johnson & Johnson, # 20391) or another disinfectant with approved material compatibility for this medical device.
- 3) 2 tanks, 1 for cleaning and 1 for disinfection - optional: 1 additional tank for rinsing with deionized/ tap water. (sufficient size for immersion of the submergible part of the probe at full length)
- 4) Soft, fluff free cloth or single use towel
- 5) Personal protective equipment (gloves, water repellent protective skirt, face protection mask or protective glasses see also instructions of the manufacturer for the detergent and the disinfectant)

### B). Bite block

- 1) Detergent: Cidezime<sup>®</sup>/ENZOL<sup>®</sup> (Johnson & Johnson, # 2258) or another cleaning agent with approved material compatibility for this medical device
- 2) Disinfectant: Cidex<sup>®</sup> OPA (Johnson & Johnson, # 20391) or another disinfectant with approved material compatibility for this medical device
- 3) 2 tanks, 1 for cleaning and 1 for disinfection - optional: 1 additional tank for rinsing with deionized/ tap water (sufficient size for immersion of the bite block)
- 4) Soft, fluff free cloth or single use towel
- 5) Personal protective equipment (gloves, water repellent protective skirt, face protection mask or protective glasses, see also instructions of the manufacturer for the detergent and the disinfectant)

## 5-5-1. Manual cleaning

### A). Probe

- 1) The temperature of the detergent solution should be between 15-30 °C [59-86 °F], concentration is 1.6%. Please note the minimum contact time of the detergent in the manufacturer's instruction. If a differing detergent is used, please also consider the approved material compatibility for this probe.
- 2) Immerse the submergible part of the probe (see figure) without connector into the detergent.
- 3) Wipe the submergible part of the probe under the surface of the detergent solution with a single-use, fluff free soft cloth to remove all visible soil. Be sure that all grooves of the probe are implemented during the cleaning process. If necessary use an appropriate cleaning brush for this purpose.
- 4) Wipe the non-submergible parts of the probe with a soft cloth dipped with a detergent. Particular attention is required for the cleaning of the lock lever and adjusting knob of the probe. Every column of the adjusting part should be wiped with a soft cloth dipped with detergent solution to assure that no soil is left in the cavities of the handling part.
- 5) Rinse the submergible part of the probe with running tap water for 1 minute.
- 6) Alternatively to step 5 suspend the submergible part of the probe in a tray filled with deionized water/tap water for 5 min.
- 7) Visually check the outer surface of the probe for cleanness. If necessary, use magnifying glass for visually check. If there is still soil visible, repeat all above steps.

### B). Bite block

- 1) The temperature of the detergent solution should be between 15-30 °C, concentration is 1.6%. Please note the minimum contact time of the detergent in the manufacturer's instruction. If a differing detergent is used, please also note the approved material compatibility for the medical device.
- 2) Immerse the bite block into the detergent.
- 3) Wipe the bite block under the surface of the detergent solution with a single-use fluff free soft cloth to remove all visible soil. Be sure that all grooves of the bite block are implemented during the cleaning process. If necessary use an appropriate cleaning brush for this purpose.
- 4) Rinse the bite block with running water 1 minute.  
(Alternatively: immerse the bite block in a tray filled with deionized/tap water for 5 minutes).
- 5) Visually check the outer surface of the bite block for cleanness. If necessary, use magnifying glass for visually check. If there is still soil visible, repeat all above steps.

## 5-5-2. Manual disinfection

## A). Probe

- 1) Before immersing the equipment, it is recommended to test the concentration of disinfectant solution before each usage. The solution Cidex® OPA is ready for use and does not need to be diluted. Test strips to verify that the appropriate concentration of Cidex® OPA is correct are available by manufacturer. Test strips will indicate a concentration above the Minimum Effective Concentration (MEC). Temperature of disinfectant solution should be minimum 20 °C[68 °F]. The minimum contact time is 5 minutes. If a differing disinfectant is used follow the manufacturer's instructions. Please also consider the material compatibility for the medical device.
- 2) Wipe the non-submergible parts of the probe with a soft and fluff free cloth with disinfectant. Particular attention is required for the disinfection of the lock lever and adjusting knob of the probe. Every column of the adjusting part should be wiped with a soft cloth dipped with disinfectant solution to assure that every cavity of the handling part is implemented.
- 3) Immerge the submergible part of the probe (see figure) into the disinfectant. Set a clock to insure the recommended contact time is observed.
- 4) Rinse the submergible part of the probe with running deionized water for 1 minute.  
(Alternatively to step 4 suspend the submergible part of the probe in a tray filled with deionized water for 5 min.)
- 5) Visually check the outer surface of the probe for that there are no leavings of the disinfectant. If necessary, repeat the rinsing.

 <b>Caution</b>
 Do not wipe the ultrasonic radiation part with alcohol. Alcohol could damage the part.

## B). Bite block

- 1) It is recommended to test the concentration of disinfectant solution before each usage. The solution Cidex® OPA is ready for use and does not need to be diluted. Test strips to verify that the appropriate concentration of Cidex® OPA is correct are available by manufacturer. Test strips will indicate a concentration above the Minimum Effective Concentration (MEC). Please also note the expiration date of the test stripes. Temperature of disinfectant solution should be minimum 20°C. The minimum contact time is 5 minutes. If a different disinfectant is used, follow the manufacturer's instructions. Please also consider the material compatibility for the medical device.
- 2) Immerge the bite block into the disinfectant. Set a clock to insure the recommended contact time is observed.
- 3) Rinse the bite block with running deionized water for 1 minute.  
(Alternatively: immerse the bite block in a tray filled with deionized water for 5 minutes).
- 4) Visually check the outer surface of the bite block for that there are no leavings of the disinfectant. If necessary, repeat the rinsing.

 <b>Warning</b>
 After finishing soaking the probe in the cleaning agent or disinfectant, thoroughly rinse it with running water (after cleaning) and deionized water (after disinfection). Residual agent can cause an adverse reaction to the operator or patient.

## 5-5-3. Cable and connector

Wipe the cable in 20 cm intervals with gauze dipped in ethyl alcohol or water, and dry it after wiping.  
Clean the connector with gauze dipped in ethyl alcohol, and dry it after cleaning.  
Clean the other parts of the probe which must not be soaked in liquid in the same manner as the connector.

 <b>Note</b>
-----------------------------------------------------------------------------------------------

If the entire length of the cable is wiped at once, a part of the cable may be wrinkled. If this occurs, pull the wrinkled part in the opposite direction to smooth it.
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## 5-6. Automated cleaning and disinfecting

## A). Probe

 <b>Warning</b>
--------------------------------------------------------------------------------------------------

 The probe cannot withstand Automated cleaning and disinfecting.
---------------------------------------------------------------------------------------------------------------------------------------------------

## B). Bite block

The following items must be provided prior to automated cleaning and disinfection.

- a) Washer disinfector according to ISO 15883 with chemo-thermal program (temperature: max 60°C [140°F])
  - b) Detergent: Korsolex® Endo-Cleaner (BODE Chemie, # 972 020) or another cleaning agent with approved material compatibility for this medical device
  - c) Disinfectant: Korsolex® Endo-Disinfectant (BODE Chemie, # 972 030) or another disinfectant with approved material compatibility for this medical device
  - d) Washer disinfector accessories: basket with lid for holding the bite block
- 1) The parameter for cleaning and disinfecting the medical device are as follows:

Program step	Water (40l)	Dosage (ml/l)	Temp. (°C)/(°F)	time (min)
Pre-Rinse	Cold water			5
Cleaning	Deionized water	5 (0.5%)	50/122	10
Rinse	Deionized water			1
Disinfection	Deionized water	10 (1%)	55/131	5
Rinse	Deionized water			1
Rinse	Deionized water		55/131	1
Drying			55/131	15

- 2) Place the bite block into the basket with lid.
- 3) Close the door of the washer disinfector and start the chemo-thermal program.
- 4) Open the door after the end of the program.
- 5) Take the bite block out of the washer disinfector and check whether it is dry. If not, proceed as described under drying.

## 5-7. Applicable cleaners and disinfectants / Suppliers List

The applicable chemical solutions are listed below.

General name	Trade name	Manufacturer
Enzyme cleaning agent	ENZOL <sup>®</sup> /Cidezime <sup>®</sup> Practical liquid 0.8V/V%	ADVANCED STERILIZATION PRODUCTS <sup>®</sup> A Johnson & Johnson company Division of Ethicon, Inc.
Alkylpolyalkylenglykolether	Korsolex <sup>®</sup> Endo-Cleaner	BODE Chemie GmbH

General name	Trade name	Manufacturer
Ortho-phthalaldehyde	CIDEX <sup>®</sup> OPA Solution 0.55%	ADVANCED STERILIZATION PRODUCTS <sup>®</sup> A Johnson & Johnson company Division of Ethicon, Inc.
Glutaral	Cidex plus <sup>®</sup>	
Glutaral	Korsolex extra *	BODE Chemie GmbH
Glutaral	Korsolex Endo- Disinfectant	BODE Chemie GmbH

Note: \* indicates that the marked disinfectant is not applicable in Canada.

## High-level disinfection

General name	Trade name	Manufacturer
Glutaraldehyde	WAVICIDE®-01 * Solution 2.65%	Medical Chemical Corporation
Glutaral	Cidex plus® Solution 3.4%	ADVANCED STERILIZATION PRODUCTS® A Johnson & Johnson company Division of Ethicon, Inc.

Note: \* indicates that the marked disinfectant is not applicable in Canada.

 <b>Warning</b>
--------------------------------------------------------------------------------------------------

- |                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  <ul style="list-style-type: none"> <li>After disinfection, thoroughly rinse the probe with deionized water.</li> <li>Residual disinfectant can cause an adverse reaction to the operator or patient.</li> </ul> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## 5-8. Drying

### A). Probe

- 1) Wipe the probe with single use, fluff free wipe or towel for removing moisture on the surface of the equipment.
- 2) If using drying heater for medical equipment, the temperature limit is a maximum of 60 °C [140 °F]. Dry until no visible moisture is left. Dry the cavities and columns of the lock lever and the adjusting knob of the handling part carefully
- 3) If using natural drying, temperature range should be between 15-30°C[59-86°F] for a minimum time of 4 hours.

### B). Bite block

- 1) Wipe the bite block with a single-use, fluff-free wipe or towel for removing moisture on the surface of the bite block.
- 2) If using a drying heater for medical equipment, the temperature limit is a maximum of 60°C [140°F]. Dry until no visible moisture is left.
- 3) If using natural drying, temperature range should be between 15-30 °C [59-86°F] for a minimum time of 4 hours.

## 5-9. Maintenance, inspection and testing

Confirm following items

- 1) the function of mechanical moving parts
- 2) the image performance when the probe is connected to the scanner
- 3) there are no abnormal exterior damages such as cracks on the surface of the equipment
- 4) Safety tests (See section 8-1)

## 5-10. Packaging

Store the disinfected probe in a dustproof environment until next application. Before sterilization it is necessary to pack all parts in a pouch suitable for sterilization, or in a tray with wrap according to ISO 11607-1 and ISO 11607-2 “Packaging for terminally sterilized devices” and ISO/TS 16775 “Packaging for terminally sterilized medical devices - Guidance on the application of ISO 11607-1 and ISO 11607-2” or the local hospital procedure. Follow the pouch manufacturer’s specifications or the local regulations for how to pack and seal the pouches. Check the sealing seam after heat sealing for any defects. In case of processing mistakes or defects the package has to be opened again and the device has to be packed and sealed again.

5-11. Sterilization

See “Table 1. Applicable cleaning, disinfection and sterilization methods” for available sterilization methods. Follow the instructions of the sterilizer manufacturer regarding usage, temperature and sterilization-time etc. Handling and maximum input to chamber of sterilizer should be according to operation manual of the sterilizer.

5-11-1. Ethylene oxide (EtO) gas sterilization

A). Probe

 <b>Warning</b>
 The probe cannot withstand ethylene oxide (EtO) gas sterilization.

B). Bite block

Sterile conditions of applicable sterilization methods are as follows.

Regarding the operation of the sterilizer, refer to the documentation supplied with the sterilizer.

Perform sterilization in the following conditions:	
Gas Type:	10% EO/ 90% HCFC
Temperature:	50 - 60°C
	122 - 140°F
Exposure Time:	More than 120 minutes
Pressurization:	162 - 200kPa
Depressurization:	13 - 8kPa
Relative humidity:	40 - 90%
Aeration is minimum	12 hours

 <b>Warning</b>
 Perform full aeration after gas sterilization. Residual gas can cause an adverse reaction to the operator or patient.

5-11-2. STERRAD® sterilization

A). Probe

⚠ Warning
⊘ The probe cannot withstand STERRAD® sterilization.

B). Bite block

Sterile conditions of applicable sterilization methods are as follows.

The applicable gas is listed below.

General name	Trade name	Manufacturer
Hydrogen peroxide (58% density)	STERRAD® Sterilization system (STERRAD® 50, 100S, 200, NX or 100NX )	ADVANCED STERILIZATION PRODUCTS® A Johnson & Johnson company Division of Ethicon, Inc.

Regarding the operation of the sterilizer, refer to the documentation supplied with the sterilizer.

Perform sterilization in the following conditions:	
STERRAD® 50, 100S or 200:	Short Cycle
Sterrad® NX or 100NX :	Standard cycle

Remark:

Some discoloration of the bite block may occur, but this does not affect performance or safety.

⚠ Caution
<p>⊘ Do not sterilize the probe using the STERRAD system if the probe is not compatible with the STERRAD system. STERRAD compatibility is shown by the STERRAD label on the connector. Perform STERRAD sterilization only for STERRAD compatible probes, otherwise it can cause damage or deterioration to the probe.</p> <p style="text-align: center;"> * STERRAD label</p>
<p>⊘ Do not put the probe directly into the sterilization pouch*. Otherwise the pouch sticks to the cable and results in damage to the cable. Completely wrap the entire probe (including the probe tip, cable and connector) with sterilization wraps* before putting it into the sterilization pouch*. *: A Johnson &amp; Johnson company Division of Ethicon, Inc. product</p>

## 5-11-3. Liquid sterilization (USA only)

## A). Probe

 <b>Warning</b>
 The probe cannot withstand liquid sterilization.

## B). Bite block

- Applicable chemical solution for sterilization  
The applicable sterilants are listed below.

General name	Trade name	Manufacturer
Glutaraldehyde	WAVICIDE®-01 * Solution 2.65%	Medical Chemical Corporation
Glutaral	Cidex plus® Solution 3.4%	ADVANCED STERILIZATION PRODUCTS® A Johnson & Johnson company Division of Ethicon, Inc.

Note: \* indicates that the marked sterilant is not applicable in Canada.

 <b>Warning</b>
 After chemical sterilization, thoroughly rinse the probe with sterile water. Residual sterilant can cause an adverse reaction to the operator or patient.

## 5-11-4. Autoclave sterilization

 <b>Warning</b>
 Please refer to " Table 1. Applicable cleaning, disinfection and sterilization methods"

## A). Probe

 <b>Warning</b>
 The probe cannot withstand autoclave sterilization.

## B). Bite block

Sterile conditions of applicable sterilization methods are as follows.  
Regarding the operation of the sterilizer, refer to the documentation supplied with the sterilizer.

Sterilize in the following conditions:    Temperature:    134°C or less

 <b>Caution</b>
 Do not carry out autoclave sterilization in a temperature condition over 134°C.

### 5-11-5. STERIS® sterilization

#### A). Probe

 <b>Warning</b>
 The probe cannot withstand STERIS® sterilization.

#### B). Bite block

The applicable product is listed below.

General name	Trade name	Manufacturer
Peracetic acid	STERIS SYSTEM 1E®	STERIS®

Regarding the operation of the sterilizer, refer to the documentation supplied with the sterilizer.

### 5-12. Storage

For details about the storage environmental conditions, see section 2-5-2 “Storage environmental conditions”.

---

## 6. Storage

---

### 6-1. Actions before storing the probe

When the equipment will not be used for an extended period of time, perform the procedures described in section 5 "Cleaning, disinfection and sterilization" and then store it in its storage case.

 <b>Caution</b>
--------------------------------------------------------------------------------------------------



The probe should be stored with the lock lever and lock knob disengaged.
--------------------------------------------------------------------------

If the probe is stored with the lock lever in the ratchet-locking position, the lever will gradually lose their ability to lock properly. This will shorten the usable life of the probe.
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### 6-2. Environmental conditions for storage

For details about the storage environmental conditions, see section 2-5-2 "Storage environmental conditions".



## 7. Moving and Transporting

### 7-1. Moving and transporting

In this section, *moving* refers to "carrying of the equipment within a facility" and *transporting* refers to "transferring using a vehicle or sending the equipment for repairs".

### 7-2. Preparing the probe and accessories for moving

Store in the storage case after performing the procedure in section 5 "Cleaning, disinfection and sterilization".

### 7-3. Packing for transportation

Store in the storage case after performing the procedure in section 5 "Cleaning, disinfection and sterilization" and then put the storage case in a cardboard box for additional protection.

### 7-4. Environmental conditions during transportation

Ambient temperature:  $-10^{\circ}\text{C}$  to  $50^{\circ}\text{C}$   
 $14^{\circ}\text{F}$  to  $122^{\circ}\text{F}$   
Relative humidity: 10% to 90%  
Atmospheric pressure: 700 hPa to 1060 hPa

#### Note

The probe is a precision equipment and is vulnerable to physical impact. Protect it by packing it properly for transportation.

Contact one of our offices and/or distributor's offices listed on the back cover when transporting the probe.



## 8. Periodic Inspection

### 8-1. Safety tests

The safety tests should be conducted at least once a year by a qualified technician. The test record should be stored for future reference.

#### Remarks 1

*Qualified technician:* personnel for conducting safety tests of medical electrical equipment.

If the user requires an appropriate qualified technician, service personnel trained by us can conduct a test at the user's expense. Contact one of our office written on the back cover.

#### Remarks 2

Make a copy of the Safety Inspection Data Sheet provided in the instruction manual of the ultrasound diagnostic instrument. Use the sheet as a test record.

### Procedure for periodic safety tests and judgment

#### (1) Test of patient leakage current from the patient connection to earth

Using the measuring instruments which usable to the requirement of IEC 60601-1 :2005, conduct the test as shown in Fig. 15 of IEC 60601-1 :2005.

Soak the probe tip in saline solution and measure the leakage current between the applied part and earth.

Do not soak probes in saline solution beyond the "IPX7 range" provided in section 2-2.

#### (2) Test of patient leakage current caused by an external voltage on the patient connection of an F-type applied part.

Using the measuring instruments which usable to the requirement of IEC 60601-1 :2005, conduct the test as shown in Fig. 16 of IEC 60601-1 :2005.

Soak the probe tip in saline solution and measure the leakage current between the applied part and earth.

Do not soak probes in saline solution beyond the "IPX7 range" provided in section 2-2.

Item	Normal condition	Single fault condition
(1) Patient leakage current from the patient connection to earth		
DC	10 $\mu$ A or less	50 $\mu$ A or less
AC	100 $\mu$ A or less	500 $\mu$ A or less
(2) Patient leakage current caused by an external voltage on the patient connection of an F-type applied part	/	5000 $\mu$ A or less

Table. Standard Values for Periodic Safety tests  
(Extract from IEC 60601-1 :2005)

### Warning



Perform a safety tests at least once a year and keep a record of the inspection results. Failure to notice an abnormal condition while using the probe can result in injury to the operator or patient. If an inspection finds an abnormal condition in the probe, immediately stop use and contact one of our offices and/or distributor's offices listed on the back cover.

## 8-2. Testing of measurement tolerances

Perform the measurements specified below using an ultrasonic phantom\* at least once per year. The test record should be stored for future reference.

- Sensitivity
- Resolution

### Remarks

Make a copy of the Measurement accuracy inspection data sheet provided in the instruction manual for the ultrasound diagnostic instrument. Use the sheet as a test record.

- \* The ultrasonic phantom is made of a substance which is similar to human tissue in terms of its response to ultrasonic waves. Regions with different textures and targets spaced at preset intervals are embedded in the phantom. Some phantoms contain a mechanism for Doppler measurement. The phantom is used to check the performance of the probe and ultrasonic diagnostic instrument, as well as to adjust the image settings.

### 8-2-1. Conducting tests

Some types of ultrasonic phantoms have targets with narrow gaps between them for confirming the resolution.

This enables you to check the level of detail that images can be viewed on the display. For phantoms with no targets, the resolution determines the fineness of the displayed textures. The sensitivity can be determined by examining the luminance of ultrasonic images. Other factors that affect the resolution include the type of connected probe, gain, focus and recording instrument. The specific testing conditions must be recorded in detail to enable proper comparison at the next inspection.

### 8-2-2. Result judgment

Compare the currently-obtained value with the value recorded at the last test. If there is a significant difference between the two values, the current value is considered to be abnormal.

It is important to note that the resolution varies depending on the type of ultrasonic phantom and phantoms generally deteriorate over time.

#### Caution



Do not use a probe or ultrasound diagnostic instrument where a problem has been found. This can result in an incorrect diagnosis. Contact one of our offices and/or distributor's offices listed on the back cover.

---

## 9. Configuration

---

### 9-1. Standard configuration

Probe	UST-52126 .....	1 set
Bite block	SP-7901 .....	2
Storage case	CB-UST9 .....	1 set
Instruction manual	MN1-5791 .....	1 copy

### 9-2. Transducer cover

CIVCO Transducer cover                      610-933

#### Remarks

The transducer cover is not included in this probe-kit.



## 10. Disposal of the Device

Recycle or dispose this equipment properly in compliance with the Waste Management and Public Cleansing Law.

### Caution



Before disposing the equipment, disinfect or take other infection-prevention measures.  
Disposal of the equipment without taking the proper preventative measures can lead to infection.

## Waste Electrical and Electronic Equipment (WEEE) Directive

The illustration on the right is required by the EU WEEE Directive to appear on all electrical and electronic equipment.

For proper disposal of this product in an EU nation, contact an EU office or agency and observe appropriate local and national regulations and laws.



■ Manufacturer



Hitachi, Ltd.

2-16-1, Higashi-Ueno, Taito-ku, Tokyo, 110-0015, Japan

■ Contact

+81-3-6284-3668

<http://www.hitachi.com/businesses/healthcare/index.html>

Overseas Offices:



Hitachi Medical Systems GmbH

Otto-von-Guericke-Ring 3 D-65205 Wiesbaden, Germany

EU Importer: Hitachi Medical Systems Europe Holding AG

Address: Sumpfstrasse 13 CH-6300 Zug, Switzerland

US Importer: Hitachi Healthcare Americas Corporation

Address: 1959 Summit Commerce Park, Twinsburg, Ohio 44087

Distributor