

# Zeus

## Diagnostic ultrasound system For General practitioner ultrasound

The MX7 ultrasound system is an essential tool for General Practitioners (GP), featuring a user-friendly interface and easy-to-use features for efficient scans. Its versatile capabilities cover abdominal, vascular, small parts, and musculoskeletal ultrasound, increasing diagnostic confidence. Seamless integration with EHR systems ensures efficient workflow and access to ultrasound images and patient records comprehensively. The MX7's lightweight design enhances portability, making it ideal for home visits.



### Simplifying ultrasound with ZST+

Be Focused on your patient and not where the focus is. Zeus has a pixel by pixel focus, where always the whole image is automatically focused.

### Ease of use with intuitive user interface iTouch

In a preset the most common settings are saved. However patients can differ and therefore settings need to be adjusted. iTouch automatically adjusts to the patient conditions very easily by just one keystroke



### Focused on diagnostic imaging

Fullscreen Mode iZoom let you concentrate on the picture and not on the information around.

### Learning tutorial on board iScan Helper

Zeus is featured with a tutorial on the system that support you to learn and train ultrasound. With dual screen function you can live refer to a reference image and elevate your ultrasound skills.

### Easily save images into your EHR

Saving images easily with network storage to a computer / server location or with DICOM to an Electronic Health Record (EHR). The data from the ultrasound system can be smoothly integrated in many environments.



Lighten  
Up infinite  
possibilities

Transducer solutions for GP



3-13 MHz



1-6 MHz

Liver in B-Mode with ZST+ pixel by pixel focused



Renal perfusion



Aorta with AMS in color mode



Thyroid in B-Mode



Ergonomic features

- Preventive magnetic power plug
- Comfortable Centric user interface layout
- Expanding probe extender
- Height adjustable cart with optional triple transducer connector available (MT2 / MT3)

mindray



M-Secure



Preventive anti-collision design

Convenient probe and cable management

Comfortable storage solutions

Slim "X"-ultra small chassis



Vibration reducing caster